


Export Bioenergy Feedstock Library Data based on lists of barcode

7/7/2020

1. From home page go to sample database under the 'Sample Database' Menu header




About Us

The Library is sponsored by DOE and hosted at Idaho National Laboratory.

<https://bioenergylibrary.inl.gov/Sample/SampleSearchGlobal.aspx>

2. Select 'GUID List Search' tab from the search menu



3. Copy and Paste GUIDs into space and select 'Search'

Home About Biomass Material Sample Database Bi

Values Search GUID List Search

59b5db91-366b-ca47-8bb5-b9961d65c4cb
ad307a3e-c068-7846-b7c2-282a46451f89
a2c183b9-d5ee-894b-a2a8-19816f984065
83753973-9a59-fb4a-9794-4496821394f3

Search Clear List

4. Check box at the top of the search results table and select 'All Samples'.

The screenshot displays a web interface for search results. At the top, there are 'Search' and 'Clear List' buttons. Below them is a section for 'Perform Group Operations on Selected Samples:' with various checkboxes for export options like 'Sample Information', 'Analyses', 'Sample Files', 'Project Files', 'Data Sets', 'Ancestry', 'Summary', 'Analysis Files', 'Citations', and 'Data Set Files'. An 'Export' button is also present. A dialog box is open in the center, asking 'Do you want to use all results, or only the current page?' with two options: 'All Samples' and 'Current Page'. The 'All Samples' option is highlighted with a red box. In the background, a table of search results is visible, with a blue checkmark in the first column of the first row highlighted by a red box. The table has columns for 'View/Edit', a green plant icon, and several columns of data including identifiers and numerical values.

View/Edit	Plant Icon	Identifier	Value 1	Value 2	Value 3
<input checked="" type="checkbox"/>		b91-366b-ca47-b9961d65c4cb	0	0	0
<input checked="" type="checkbox"/>		ad307a3e-c068-7846-b7c2-282a46451f89	0	0	0
<input checked="" type="checkbox"/>		a2c183b9-d5ee-894b-a2a8-19816f984065	0	0	0
<input checked="" type="checkbox"/>		83753973-9a59-fb4a-9794-4496821394f3	0	0	0

5. Check boxes for export options and select 'Export'

Search Clear List

Perform Group Operations on Selected Samples:

Sample Information Analyses Sample Files Project Files Data Sets
Export Ancestry Summary Analysis Files Citations Data Set Files **Export**

Drag a column header and drop it here to group by that column

Page size: 100

<input checked="" type="checkbox"/>			Analysis	Hierarchy	Files	Citations	Data Sets	GUID Id
<input checked="" type="checkbox"/>		View/Edit			0	0	0	59b5db9 8bb5-b9
<input checked="" type="checkbox"/>		View/Edit			0	0	0	ad307a1 b7c2-28
<input checked="" type="checkbox"/>		View/Edit			0	0	0	a2c183t a2a8-19
<input checked="" type="checkbox"/>		View/Edit			0	0	0	837539t 9794-44

Page size: 100

*Note that selecting 'Sample Files', 'Project Files', 'Analysis Files', or 'Data Set Files' will result in the export of a Zip File with the selected export option associated files.